

Anti-T3 antibody [10550] (STJ400267)

STJ400267

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-T3 is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	10550
Concentration	> 1.0 mg/ml
Conjugation	Unconjugated
Purification	Greater than 95% purity
Dilution Range	
Formulation	0.1 M PBS, pH 7.4, 0.9 % NaCl, 0.05 % NaN ₃ as a preservative
Isotype	IgG1
Storage	60 months from manufacturing at 2–8°C. Temperature analysis:-70 °C, 21 days Not Determined (N/D)-20 °C, 14 days OK +4 °C, 14 days OK +35 °C, 14 days OK +45 °C, 7 days OK. Shelf life stability testing at 2–8 °C in the product buffer.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human triiodothyronine (T3)
Immunogen Region	
Specificity	
Immunogen Sequence	