

Anti-Vitamin D antibody [3850] (STJ400268)

STJ400268

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Vitamin D is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Vitamin D

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3850
Concentration	> 1.0 mg/ml
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.1 M PBS, pH 7.4, 0.9 % NaCl, 0.05 % NaN ₃ as a preservative
Isotype	IgG1
Storage Instruction	60 months from manufacturing at 2–8°C.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Natural 25 OH Vitamin D ₂ and 25 OH Vitamin D ₃
Immunogen Region	
Specificity	
Immunogen Sequence	