

## Anti-Calcitonin antibody [4050] (STJ400270)

STJ400270

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Calcitonin is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	4050
<b>Concentration</b>	> 1.0 mg/ml
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Greater than 95% purity
<b>Dilution Range</b>	
<b>Formulation</b>	0.1 M BBS, pH 8, 0.9 % NaCl
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	60 months from manufacturing at 2–8°C.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Human calcitonin
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	