

Anti-PTH antibody (57-61) [10750] (STJ400272)

STJ400272

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-PTH (57-61) is suitable for use in ELISA research applications. |
| Applications | ELISA |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | 10750 |
| Concentration | > 1.0 mg/ml |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | 0.1 M PBS, pH 7.4, 0.9 % NaCl, 0.05 % NaN ₃ as a preservative |
| Isotype | IgG1 |
| Storage Instruction | 60 months from manufacturing at 2–8°C. |

TARGET INFORMATION

| | |
|---------------------------|---------------------------------|
| Gene ID | 138428 |
| Gene Symbol | PTRH1 |
| Uniprot ID | PTH_HUMAN |
| Immunogen | Human parathyroid hormone (PTH) |
| Immunogen Region | 57-61 |
| Specificity | |
| Immunogen Sequence | |