

**Anti-NT-proBNP antibody [RC1306] (STJ400297)**

STJ400297

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse, Human monoclonal antibody anti-NT-proBNP is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Mouse, Human
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	RC1306
<b>Concentration</b>	5.0 mg/ml (+/-10 %) (A280 nm, 1 mg/ml, 1 cm = 1.6)
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
<b>Isotype</b>	human IgG1
<b>Storage Instruction</b>	Unspecified, storage at 2–8 °C

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	Antibody recognizes human N-terminal pro-B-type natriuretic peptide
<b>Immunogen Sequence</b>	