

Anti-NT-proBNP antibody [RC1307] (STJ400298)

STJ400298

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse, Human monoclonal antibody anti-NT-proBNP is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse, Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RC1307
Concentration	5.0 mg/ml (+/-10 %) , (A280nm, 1 mg/ml, 1 cm= 1.4)
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % Na ₃ as a preservative
Isotype	human IgG1
Storage Instruction	Unspecified, storage at 2–8 °C

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	Antibody recognizes human N-terminal pro-B-type natriuretic peptide
Immunogen Sequence	