

Anti-Prealbumin antibody [11604] (STJ400302)

STJ400302

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-Prealbumin is suitable for use in ELISA research applications. |
| Applications | ELISA |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | 11604 |
| Concentration | 5.0 mg/ml (+/-10 %) |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative |
| Isotype | IgG1 |
| Storage Instruction | Unspecified, storage at 2-8 °C |

TARGET INFORMATION

| | |
|---------------------------|--------------------------------------|
| Gene ID | 7276 |
| Gene Symbol | TTR |
| Uniprot ID | TTHY_HUMAN |
| Immunogen | |
| Immunogen Region | |
| Specificity | Antibody recognizes human prealbumin |
| Immunogen Sequence | |