

Anti-SAA antibody [2212] (STJ400312)

STJ400312

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-SAA is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2212
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % Na ₃ N as a preservative
Isotype	IgG1
Storage Instruction	Unspecified, storage at 2–8 °C

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	Antibody recognizes human serum amyloid A (SAA)
Immunogen Sequence	