

Anti-SHBG antibody [6007] (STJ400318)

STJ400318

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-SHBG is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	6007
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.9 % NaCl, 0.095 % NaN ₃ as a preservative
Isotype	IgG1
Storage Instruction	24 months from manufacturing at 2–8 °C

TARGET INFORMATION

Gene ID	6462
Gene Symbol	SHBG
Uniprot ID	SHBG_HUMAN
Immunogen	
Immunogen Region	
Specificity	Antibody recognizes human sex hormone-binding globulin
Immunogen Sequence	