

Anti-SARS-CoV-2 N antibody [3510] (STJ400323)

STJ400323

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-SARS-CoV-2 N is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3510
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Unspecified, storage at 2-8 °C

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	Antibody recognizes SARS-CoV-2 nucleocapsid protein
Immunogen Sequence	