

Anti-Thyroxine antibody [6902] (STJ400326)

STJ400326

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Thyroxine is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Thyroxine

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	6902
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
Isotype	IgG1
Storage Instruction	18 months from manufacturing at 2–8 °C

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	Antibody recognizes human thyroxine
Immunogen Sequence	