

Anti-THC antibody (STJ400332)

STJ400332

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Cannabinoid is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Cannabinoid

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	≥ 1.5 mg/ml
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.05M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl, 0.05% NaN ₃ as a preservative
Isotype	IgG1
Storage Instruction	12 months from manufacturing at 2–8 °C

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	Antibody recognizes THC and its metabolites
Immunogen Sequence	