

Anti-E3 antibody [RC10601] (STJ400338)

STJ400338

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse, Human monoclonal antibody anti-Estriol is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse, Human
Reactivity	Estriol

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RC10601
Concentration	5.0 mg/ml (+/-10 %) (A280nm, 1 mg/ml, 1 cm = 1.5)
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	human IgG1
Storage Instruction	Unspecified, storage at 2–8 °C

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	Antibody recognizes estriol
Immunogen Sequence	