

Anti-HTR2A antibody (400-450) {Biotin} (STJ500007)

STJ500007

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-HTR2A (400-450) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:100 ICC: 1:50-1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3356
Gene Symbol	HTR2A
Uniprot ID	5HT2A_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-450 on human 5-hydroxytryptamine receptor 2A.
Immunogen Region	400-450
Specificity	
Immunogen Sequence	