

## Anti-TNFAIP3 antibody (690-790) (STJ500012)

STJ500012

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-TNFAIP3 (690-790) is suitable for use in ELISA and Western Blot research applications.
<b>Applications</b>	ELISA, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Monkey, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.56-0.62 µg/µl
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 ELISA: 1:10, 000 IP: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	7128
<b>Gene Symbol</b>	TNFAIP3
<b>Uniprot ID</b>	TNAP3_HUMAN
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 690-790 on human Tumor necrosis factor alpha-induced protein 3.
<b>Immunogen Region</b>	690-790
<b>Specificity</b>	
<b>Immunogen Sequence</b>	



Western blot with Anti-TNFAIP3 Antibody (STJ500012) with peptide on at 1:500 dilution in buffer.