

Anti-ABHD2 antibody (200-250) (STJ500015)

STJ500015

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ABHD2 (200-250) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54608
Gene Symbol	Abhd2
Uniprot ID	ABHD2_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 200-250 on mouse a/b hydrolase II protein.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	



Immunohistochemistry of mouse kidney with Anti-ABHD2 primary antibody (STJ500015) at 1:50 dilution in Buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.