

Anti-ABRAXAS1 antibody (350-400) (STJ500021)

STJ500021

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ABRAXAS1 (350-400) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.52-0.63 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84142
Gene Symbol	ABRAXAS1
Uniprot ID	ABRX1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 350-400 on human BRCA1-A complex subunit Abraxas protein.
Immunogen Region	350-400
Specificity	
Immunogen Sequence	



Western blot with Anti-ABRAXAS1 Antibody (STJ500021) with peptide on at 1:500 dilution in buffer.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081