

Anti-TNK2 antibody (250-300) (STJ500027)

STJ500027

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TNK2 (250-300) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.59-0.68 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10188
Gene Symbol	TNK2
Uniprot ID	ACK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-400 on ACK1 protein.
Immunogen Region	250-300
Specificity	
Immunogen Sequence	



Western blot Immunoprecipitation with Anti-TNK2 Antibody (STJ500027) of solubilized COS1 cells expressing constitutively active ACK1, at 1:500 dilution in buffer.