

Anti-ACKR2 antibody (200-300) (STJ500030)

STJ500030

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ACKR2 (200-300) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1238
Gene Symbol	ACKR2
Uniprot ID	ACKR2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 200-300 on ACRD6 protein.
Immunogen Region	200-300
Specificity	
Immunogen Sequence	



Western blot with Anti-ACKR2 Antibody (STJ500030) of rat lung at 1:500 dilution in buffer.