

Anti-ACTN4 antibody (140-190) (STJ500045)

STJ500045

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ACTN4 (140-190) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.62 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	81
Gene Symbol	ACTN4
Uniprot ID	ACTN4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 140-190 on human alpha Actinin 4 protein.
Immunogen Region	140-190
Specificity	
Immunogen Sequence	



Western blot with Anti-ACTN4 Antibody (STJ500045) of
1) rat podocyte membranes and 2) NIH3T3 cells.
Antibody was diluted 1:500 in buffer.