

## Anti-AFP antibody (Full Length) {Biotin} (STJ500052)

STJ500052

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-AFP (Full Length) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IHC, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Monkey, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.68 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 ELISA: 1:10,000 IP: 1:200 IHC: 1:100
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	174
<b>Gene Symbol</b>	AFP
<b>Uniprot ID</b>	FETA_HUMAN
<b>Immunogen</b>	Full length native protein (purified) corresponding to Human alpha 1 Fetoprotein. Purified human alpha 1 Fetoprotein.
<b>Immunogen Region</b>	Full Length
<b>Specificity</b>	
<b>Immunogen Sequence</b>	