

Anti-ALDH1A1 antibody (441-491) {Biotin} (STJ500064)

STJ500064

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-ALDH1A1 (441-491) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications. |
| Applications | ELISA, IHC, IP, WB |
| Host/Source | Rabbit |
| Reactivity | Mouse, Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.5 µg/µl |
| Conjugation | Biotin |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:100 |
| Formulation | |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 24188 |
| Gene Symbol | Aldh1a1 |
| Uniprot ID | AL1A1_RAT |
| Immunogen | Synthetic peptide taken within amino acid region 441-491 on rat retinal dehydrogenase 1 protein. |
| Immunogen Region | 441-491 |
| Specificity | |
| Immunogen Sequence | |