

Anti-ALDH1A1 antibody (441-491) {FITC} (STJ500065)

STJ500065

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ALDH1A1 (441-491) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	24188
Gene Symbol	Aldh1a1
Uniprot ID	AL1A1_RAT
Immunogen	Synthetic peptide taken within amino acid region 441-491 on rat retinal dehydrogenase 1 protein.
Immunogen Region	441-491
Specificity	
Immunogen Sequence	