

Anti-AMPK alpha 2 antibody (480-530) {Biotin} (STJ500079)

STJ500079

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-AMPK alpha 2 (480-530) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:250
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	78975
Gene Symbol	Prkaa2
Uniprot ID	AAPK2_RAT
Immunogen	Synthetic peptide taken within amino acid region 480-530 on rat 5 NA-AMP-activated protein kinase catalytic subunit alpha-2 protein.
Immunogen Region	480-530
Specificity	
Immunogen Sequence	