

Anti-ANGPTL4 antibody (50-100) (STJ500081)

STJ500081

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ANGPTL4 (50-100) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57875
Gene Symbol	Angptl4
Uniprot ID	ANGL4_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 50-100 on mouse Angiopoietin-related protein 4.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	



Western blot with Anti-Angiopoietin-like 4 Antibody (STJ500081) of mouse at dilution in buffer. WT) Wild type KO) ANGPTL4 knockout mouse macrophages.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081