

Anti-ANGPTL5 antibody (340-388) {Biotin} (STJ500085)

STJ500085

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ANGPTL5 (340-388) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.66 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	253935
Gene Symbol	ANGPTL5
Uniprot ID	ANGL5_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 340-388 on human ANGPTL5 protein.
Immunogen Region	340-388
Specificity	
Immunogen Sequence	