

Anti-ANXA10 antibody (60-100) {Biotin} (STJ500094)

STJ500094

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ANXA10 (60-100) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11199
Gene Symbol	ANXA10
Uniprot ID	ANX10_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 60-100 on human Annexin A10 protein.
Immunogen Region	60-100
Specificity	
Immunogen Sequence	