

Anti-ANXA6 antibody (500-560) (STJ500096)

STJ500096

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ANXA6 (500-560) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	309
Gene Symbol	ANXA6
Uniprot ID	ANXA6_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 500-560 on human Annexin A6 protein.
Immunogen Region	500-560
Specificity	
Immunogen Sequence	



Western blot with Anti-ANXA6 Antibody (STJ500096) with peptide on at 1:500 dilution in buffer.