

Anti-ANXA9 antibody (1-100) (STJ500099)

STJ500099

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ANXA9 (1-100) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8416
Gene Symbol	ANXA9
Uniprot ID	ANXA9_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-100 on human Annexin A9 protein.
Immunogen Region	1-100
Specificity	
Immunogen Sequence	



Immunohistochemistry of baboon lung with Anti-ANXA9 primary antibody (STJ500099) at 1:50 dilution in Buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40x.