

Anti-APLN antibody (25-77) (STJ500102)

STJ500102

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-APLN (25-77) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8862
Gene Symbol	APLN
Uniprot ID	APEL_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 25-77 on human Apelin protein.
Immunogen Region	25-77
Specificity	
Immunogen Sequence	



Immunohistochemistry of rat small intestine with Anti-APLN primary antibody (STJ500102) at 1:50 dilution in Buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.