

**Anti-APH1B antibody (235-257) {Biotin} (STJ500106)**

STJ500106

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-APH1B (235-257) is suitable for use in ELISA and Western Blot research applications.
<b>Applications</b>	ELISA, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	C.elegans, Human, Mouse, Rat

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.6 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	83464
<b>Gene Symbol</b>	APH1B
<b>Uniprot ID</b>	APH1B_HUMAN
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 235-257 of Human Presenilin stabilizing protein.
<b>Immunogen Region</b>	235-257
<b>Specificity</b>	APH1B-BIOTIN does not cross react to other components of gamma secretase complex nor with other proteins.
<b>Immunogen Sequence</b>	