

Anti-APLNR antibody (150-200) {Biotin} (STJ500109)

STJ500109

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-APLNR (150-200) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	187
Gene Symbol	APLNR
Uniprot ID	APJ_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 150-200 human Apelin receptor protein.
Immunogen Region	150-200
Specificity	
Immunogen Sequence	