

## Anti-AQP11 antibody (140-190) {Biotin} (STJ500121)

STJ500121

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-AQP11 (140-190) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	CM, ELISA, ICC, IF, IHC, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Pig, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.57-0.67 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:500 ICC: 1:50 IF: 1:50
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	286758
<b>Gene Symbol</b>	Aqp11
<b>Uniprot ID</b>	AQP11_RAT
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 140-190 on Aquaporin 11 protein.
<b>Immunogen Region</b>	140-190
<b>Specificity</b>	
<b>Immunogen Sequence</b>	