

Anti-AQP11 antibody (140-190) {FITC} (STJ500122)

STJ500122

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-AQP11 (140-190) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Pig, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.57-0.67 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:500 ICC: 1:50 IF: 1:50
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	286758
Gene Symbol	Aqp11
Uniprot ID	AQP11_RAT
Immunogen	Synthetic peptide taken within amino acid region 140-190 on Aquaporin 11 protein.
Immunogen Region	140-190
Specificity	
Immunogen Sequence	