

Anti-AQP3 antibody (200-250) {Biotin} (STJ500130)

STJ500130

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-AQP3 (200-250) is suitable for use in ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Bovine, Human, Mouse, Rat, Sheep

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50-0.57 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100 ICC: 1:100 IF: 1:100
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	65133
Gene Symbol	Aqp3
Uniprot ID	AQP3_RAT
Immunogen	Synthetic peptide taken within amino acid region 200-250 on rat Aquaporin 3 protein.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	