

Anti-AQP6 antibody (220-270) {Biotin} (STJ500139)

STJ500139

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-AQP6 (220-270) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Bovine, Human, Mouse, Rat, Sheep

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65-0.75 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:400 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29170
Gene Symbol	Aqp6
Uniprot ID	AQP6_RAT
Immunogen	Synthetic peptide taken within amino acid region 220-270 on rat Aquaporin 6 protein.
Immunogen Region	220-270
Specificity	
Immunogen Sequence	