

Anti-AQR antibody (1250-1300) {Biotin} (STJ500151)

STJ500151

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-AQR (1250-1300) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.57-0.60 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:250
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9716
Gene Symbol	AQR
Uniprot ID	AQR_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1250-1300 on human RNA helicase aquarius protein.
Immunogen Region	1250-1300
Specificity	
Immunogen Sequence	