

Anti-ARF1 antibody {FITC} (STJ500155)

STJ500155

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ARF1 is suitable for use in ELISA, Immunocytochemistry, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, ICC, IF, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.58-0.65 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 ICC: Use at an assay dependent concentration
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	64310
Gene Symbol	Arf1
Uniprot ID	ARF1_RAT
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on Rat ARF1 protein. The peptide was covalently modified post-synthetically to achieve desired antigenicity.
Immunogen Region	
Specificity	
Immunogen Sequence	