

Anti-ARHGEF7 antibody (595-805) {FITC} (STJ500161)

STJ500161

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ARHGEF7 (595-805) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.59 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10,000 ELISA: 1:10,000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8874
Gene Symbol	ARHGEF7
Uniprot ID	ARHG7_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 595-805 on human ARHGEF7 protein.
Immunogen Region	595-805
Specificity	
Immunogen Sequence	