

## Anti-CYP19A1 antibody (400-450) (STJ500162)

STJ500162

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-CYP19A1 (400-450) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IHC, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5-1.5 µg/µl
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:150
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	1588
<b>Gene Symbol</b>	CYP19A1
<b>Uniprot ID</b>	CP19A_HUMAN
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 400-450 on human Aromatase protein.
<b>Immunogen Region</b>	400-450
<b>Specificity</b>	
<b>Immunogen Sequence</b>	



Western blot with Anti-CYP19A1 Antibody (STJ500162) of mouse brain at 1:500 dilution in buffer.