

Anti-SLC7A10 antibody {Biotin} (STJ500167)

STJ500167

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SLC7A10 is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50-0.75 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:1,000 DB: 1:4,000 ELISA: 1:4,000 IP: 1:200 IHC: 1:50 ICC: 1:50 IF: 1:50 CFM: 1:50
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	56301
Gene Symbol	SLC7A10
Uniprot ID	AAA1_HUMAN
Immunogen	C-epitope Synthetic peptide corresponding to unique amino acid sequence on Na-independent D-amino acid transporter protein.
Region	
Specificity	
Immunogen Sequence	