

Anti-ASIC2 antibody (450-512) (STJ500174)

STJ500174

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ASIC2 (450-512) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	40
Gene Symbol	ASIC2
Uniprot ID	ASIC2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 450-512 on human Acid-sensing ion channel 2 protein.
Immunogen Region	450-512
Specificity	
Immunogen Sequence	



Western blot with Anti-ASIC2 Antibody (STJ500174) of mouse kidney at 1:1000 dilution in buffer.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081