

Anti-ASNA1 antibody (300-348) {Biotin} (STJ500178)

STJ500178

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ASNA1 (300-348) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.52-.53 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	439
Gene Symbol	GET3
Uniprot ID	GET3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-348 on human ATPase ASNA1 protein.
Immunogen Region	300-348
Specificity	
Immunogen Sequence	