

Anti-AURKB antibody (300-344) {Biotin} (STJ500193)

STJ500193

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-AURKB (300-344) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications. |
| Applications | ELISA, IP, WB |
| Host/Source | Rabbit |
| Reactivity | Human, Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.5-1.5 µg/µl |
| Conjugation | Biotin |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:250 |
| Formulation | |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 9212 |
| Gene Symbol | AURKB |
| Uniprot ID | AURKB_HUMAN |
| Immunogen | Synthetic peptide taken within amino acid region 300-344 on human Aurora Kinase B protein. |
| Immunogen Region | 300-344 |
| Specificity | |
| Immunogen Sequence | |