

Anti-AXIN2 antibody (800-843) {Biotin} (STJ500196)

STJ500196

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-AXIN2 (800-843) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8313
Gene Symbol	AXIN2
Uniprot ID	AXIN2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 800-843 on human Axin 2 protein.
Immunogen Region	800-843
Specificity	
Immunogen Sequence	