

Anti-Axl antibody (840-880) {Biotin} (STJ500199)

STJ500199

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Axl (840-880) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.59-0.65 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	26362
Gene Symbol	Axl
Uniprot ID	UFO_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 840-880 on mouse AXL protein. The peptide was post translationally modified to achieve proper orientation and antigenicity.
Immunogen Region	840-880
Specificity	
Immunogen Sequence	