

Anti-BACE2 antibody (350-396) {Biotin} (STJ500211)

STJ500211

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-BACE2 (350-396) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.59 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25825
Gene Symbol	BACE2
Uniprot ID	BACE2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 350-396 on human BACE2B protein.
Immunogen Region	350-396
Specificity	This antibody will only detect BACE2B. This antibody will not cross react with BACE2A or BACE2B.
Immunogen Sequence	