

Anti-BCKDK antibody (170-200) {Biotin} (STJ500214)

STJ500214

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-BCKDK (170-200) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10295
Gene Symbol	BCKDK
Uniprot ID	BCKD_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 170-200 on human BCKD protein.
Immunogen Region	170-200
Specificity	
Immunogen Sequence	