

Anti-BCL6 antibody (650-700) (STJ500220)

STJ500220

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-BCL6 (650-700) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.6 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	604
Gene Symbol	BCL6
Uniprot ID	BCL6_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 650-700 on human B-cell lymphoma 6 protein isoform 1.
Immunogen Region	650-700
Specificity	Will detect BCL6 isoform a, BCL6 isoform b and BCL6S protein.
Immunogen Sequence	



Immunohistochemistry of rat spleen with Anti-BCL6 primary antibody (STJ500220) at 1:50 dilution in Buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.