

Anti-BMP1 antibody (1-50) {Biotin} (STJ500266)

STJ500266

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-BMP1 (1-50) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50-0.56 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:2,000 ELISA: 1:2,000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	649
Gene Symbol	BMP1
Uniprot ID	BMP1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human BMP1 protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	